

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Agilent (HP)	1080	Bare Lamp	D001
Agilent (HP)	1081	Bare Lamp	D001
Agilent (HP)	1081B	Bare Lamp	D001
Agilent (HP)	1084B	Bare Lamp	D001
Altex / Hitachi	100 (HPLC)		D002
Altex	155		D012
Altex	163		D016
Altex	164	Bare Lamp	D011
Altex	165	Bare Lamp	D011
Altex	167	Bare Lamp	D011
Aminco	DW2C	Check Terminations	D002
Barspec	All Models		D002
Bausch & Lomb	200	Use clip from old lamp	D017
Bausch & Lomb	700	Use clip from old lamp	D017
Bausch & Lomb	710	Use clip from old lamp	D017
Bausch & Lomb	2000	Use clip from old lamp	D017
Bausch & Lomb	3100		D016
Bausch & Lomb	3200		D016
Bausch & Lomb	Aquamate UV		D006
Bausch & Lomb	DE 50	Use clip from old lamp	D017
Bausch & Lomb	Genesys 8		D006
Bausch & Lomb	Spectronic 21	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 21D	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 501	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 601	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 1001	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 1001 plus	Use clip from old lamp	D017
Bausch & Lomb	Spectronic 1201	Use clip from old lamp	D017
Beckman	20		D011
Beckman	24		D011
Beckman	25		D011
Beckman	26		D011
Beckman	34		D011
Beckman	35		D011
Beckman	36		D011
Beckman	100	Check Terminations	D002
Beckman	155		D012
Beckman	163		D016
Beckman	164		D011
Beckman	165		D011
Beckman	167		D011
Beckman	5200		D011
Beckman	Acta		D011
Beckman	DU6		D011
Beckman	DU7		D011
Beckman	DU20		D011
Beckman	DU30		D011
Beckman	DU37		D011

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Beckman	DU40		D011
Beckman	DU50		D011
Beckman	System Gold	Check Terminations	D011
Biotronic	BT3035	Bare Lamp	D012
Cecil Instruments	1000 Series	Pt. No 2202-0142	D007
Cecil Instruments	2000 Series	Pt. No 2202-0142	D007
Cecil Instruments	3000 Series	Pt. No 2202-0142	D007
Cecil Instruments	5000 Series	Pt. No 2202-0142	D007
Cecil Instruments	6000 Series	Pt. No 2202-0142	D007
Cecil Instruments	7000 Series	Pt. No 2202-0142	D007
Cecil Instruments	8000 Series	Pt. No 2202-0142	D007
Cecil Instruments	9000 Series	Pt. No 2202-0142	D007
Cecil Instruments	CD 4201	Pt. No 2202-0142	D007
Cecil Instruments	CD 4400	Pt. No 2202-0142	D007
Cecil Instruments	CD 6600S	Pt. No 2202-0142	D007
Cecil Instruments	CE 1010	Pt. No 2202-0142	D007
Cecil Instruments	CE 1020	Pt. No 2202-0142	D007
Cecil Instruments	CE 1070	Pt. No 2202-0142	D007
Cecil Instruments	CE 1200	Pt. No 2202-0142	D007
Cecil Instruments	CE 1220	Pt. No 2202-0142	D007
Cecil Instruments	CE 2021	Pt. No 2202-0142	D007
Cecil Instruments	CE 2202	Pt. No 2202-0142	D007
Cecil Instruments	CE 2272	Pt. No 2202-0142	D007
Cecil Instruments	CE 2292	Pt. No 2202-0142	D007
Cecil Instruments	CE 4200	Pt. No 2202-0142	D007
Cecil Instruments	CE 4400	Pt. No 2202-0142	D007
Cecil Instruments	CE 5501	Pt. No 2202-0142	D007
Cecil Instruments	CE 5502	Pt. No 2202-0142	D007
Cecil Instruments	CE 6000	Pt. No 2202-0142	D007
Cecil Instruments	CE 6602	Pt. No 2202-0142	D007
Cecil Instruments	CE 6603	Pt. No 2202-0142	D007
Cecil Instruments	CE 7200	Pt. No 2202-0142	D007
Cecil Instruments	CE 8020	Pt. No 2202-0142	D007
Cecil Instruments	CE 9020	Pt. No 2202-0142	D007
Cecil Instruments	CE 9050	Pt. No 2202-0142	D007
Cecil Instruments	CE 9080	Pt. No 2202-0142	D007
Cecil Instruments	CE 9090	Pt. No 2202-0142	D007
Cecil Instruments	Series 2 Models	Pt. No 2202-0142	D007
Coleman	55	Check Terminations	D002
Coleman	124	Check Terminations	D002
Coleman	550	Check Terminations	D002
Coleman	552	Check Terminations	D002
Constametric	Spectromonitor III		D016
Dupont	825		D002
Dupont	850		D002
EM Science	All	Use old connector	D002
Erma	ERC 7210	Use old connector	D002

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Femto	700		D017
Femto	700 Plus		D017
Femto	700S		D017
Femto	800XL		D017
GBC	906		D008
GBC	All AA Instruments		D008
GBC	Cintra 10		D004
GBC	Cintra 20		D004
GBC	Cintra 40		D004
GBC	Cintra 101		D004
GBC	Cintra 202		D004
GBC	Cintra 303		D004
GBC	Cintra 404		D004
GBC	UV 911		D008
GBC	UV VIS 914		D004
GBC	UV VIS 916		D004
GBC	UV VIS 918		D004
GBC	UV VIS 920		D004
Gilford	200		D001
Gilford	220		D001
Gilford	222		D001
Gilford	230		D001
Gilford	240		D001
Gilford	250		D001
Gilford	252		D001
Gilford	260		D001
Gilford	2000		D001
Gilford	2400		D001
Gilford	2600		D001
Gilford	Response		D002
Hewlett Packard	1080	Use bracket from old lamp	D001
Hewlett Packard	1081	Use bracket from old lamp	D001
Hewlett Packard	1081B	Use bracket from old lamp	D001
Hewlett Packard	1082B	Use bracket from old lamp	D001
Hewlett Packard	1084	Use bracket from old lamp	D001
Hewlett Packard	1084B	Use bracket from old lamp	D001
Hitachi	101		D012
Hitachi	102		D012
Hitachi	111		D012
Hitachi	124		D002
Hitachi	181	Prealigned	D005
Hitachi	200	Prealigned	D005
Hitachi	220	Prealigned	D005
Hitachi	220A	Prealigned	D005
Hitachi	300	Prealigned	D005
Hitachi	330	Prealigned	D005
Hitachi	330LC	Prealigned	D005
Hitachi	340	Prealigned	D005

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Hitachi	100-10		D002
Hitachi	100-30		D002
Hitachi	100-40		D002
Hitachi	100-50		D002
Hitachi	100-60		D002
Hitachi	150-20	Prealigned	D005
Hitachi	655A22	Prealigned	D005
Hitachi	L2000	Prealigned	D005
Hitachi	L2500	Prealigned	D005
Hitachi	L3000	Prealigned	D005
Hitachi	L3200	Prealigned	D005
Hitachi	L4000	Prealigned	D005
Hitachi	L4200	Prealigned	D005
Hitachi	L4225	Prealigned	D005
Hitachi	L4250	Prealigned	D005
Hitachi	L4500	Prealigned	D005
Hitachi	L4520	Prealigned	D005
Hitachi	L4720	Prealigned	D005
Hitachi	L6000	Prealigned	D005
Hitachi	L7000	Prealigned	D005
Hitachi	L7400	Prealigned	D005
Hitachi	L7420	Prealigned	D005
Hitachi	L7450	Prealigned	D005
Hitachi	L7450A	Prealigned	D005
Hitachi	L7455 DAD Lachrom	Prealigned	D005
Hitachi	Lachrom 655A detector	Prealigned	D005
Hitachi	U1000	Prealigned	D005
Hitachi	U1100	Prealigned	D005
Hitachi	U1500	Prealigned	D005
Hitachi	U2000	Prealigned	D005
Hitachi	U2001	Prealigned	D005
Hitachi	U3010	Prealigned	D005
Hitachi	U3210	Prealigned	D005
Hitachi	U3310	Prealigned	D005
Hitachi	U3410	Prealigned	D005
Hitachi	U3501	Prealigned	D005
Hitachi	U4000	Prealigned	D005
Hitachi	U4001	Prealigned	D005
Hitachi	U4010	Prealigned	D005
Hitachi	UL1100	Prealigned	D005
Hitachi / Merck	U-1800		D015
Hitachi / Merck	U-2800		D015
Hitachi / Merck	U-2810		D015
IBM Instruments	Detectors		D016
ICI Instruments	LC1200	Check Terminations	D002
Isco	V4		D001
Jasco	Uvidec 1		D012
Jasco	Uvidec 100		D012

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Jasco	Uvidec 100 II		D012
Jasco	Uvidec 2		D012
Jasco	Uvidec 320		D012
Jasco	Uvidec 340		D012
Jasco	Uvidec 4		D012
Jasco	Uvidec 420		D012
Jasco	Uvidec 4D		D002
Jasco	Uvidec 505		D012
Jasco	Uvidec 510		D012
Jasco	Uvidec 710		D012
Jasco	V530		D003
Jasco	V550		D003
Jasco	V650		D003
Jasco	V660		D003
Jasco	V670		D003
Joyce LoebL	Scan 400		D001
Kontron	Uvikon 725		D002
Kontron	725 LC		D002
LDC	3000		D016
LDC	3100		D016
LDC	3200		D016
LDC	4000 Series		D016
LDC	4100		D016
LDC	4200		D016
LDC	5000		D016
LDC	Spectromonitor 1		D016
LDC	Spectromonitor 2		D016
LDC	Spectromonitor 3		D016
LDC	Spectromonitor 3000		D016
LDC	Spectromonitor 3200		D016
LDC	Spectromonitor 4100		D016
LDC	Spectromonitor D		D016
LKB Pharmacia	4050	Resolder old circuit board	D002
LKB Pharmacia	4054	Resolder old circuit board	D002
LKB Pharmacia	TDS S/N>205	Bare Lamp	D002
LKB Pharmacia	Ultraspec II(S/N>6598)	Resolder old circuit board	D002
LKB Pharmacia	Ultraspec III		D002
LKB Pharmacia	Ultraspec Plus (S/N>175)		D002
LKB Pharmacia	Ultraspec+(S/N>175)		D002
Merck / Hitachi	655A/22	Prealigned	D005
Merck / Hitachi	L3000	Prealigned	D005
Merck / Hitachi	L4000	Prealigned	D005
Merck / Hitachi	L4500	Prealigned	D005
Merck / Hitachi	L4520	Prealigned	D005
Merck / Hitachi	L4720	Prealigned	D005
Merck / Hitachi	L7400	Prealigned	D005
Merck / Hitachi	L7420	Prealigned	D005
Merck / Hitachi	L7450	Prealigned	D005

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Merck / Hitachi	L7455 DAD	Prealigned	D005
Merck / Hitachi	U-1800	Prealigned	D015
Merck / Hitachi	U-2800	Prealigned	D015
Merck / Hitachi	U-2810	Prealigned	D015
Micromeritics	787		D009
Micromeritics	788		D009
Micronal	B582		D003
Milton Roy	3100		D016
Milton Roy	3200		D016
Milton Roy - Spectronic	200		D017
Milton Roy - Spectronic	700		D017
Milton Roy - Spectronic	710		D017
Milton Roy - Spectronic	2000		D017
Milton Roy - Spectronic	3100		D016
Milton Roy - Spectronic	3200		D016
Milton Roy - Spectronic	DE 50		D017
Milton Roy - Spectronic	Spectronic 21		D017
Milton Roy - Spectronic	Spectronic 21D		D017
Milton Roy - Spectronic	Spectronic 501		D017
Milton Roy - Spectronic	Spectronic 601		D017
Milton Roy - Spectronic	Spectronic 1001		D017
Milton Roy - Spectronic	Spectronic 1001 plus		D017
Milton Roy - Spectronic	Spectronic 1201		D017
NICOLET	Evolution 100	Prealigned	D006
NICOLET	Evolution 500	Prealigned	D006
NICOLET	Evolution 600	Prealigned	D006
Perkin Elmer / Hitachi	124	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	139	Bare lamp/check	D002
Perkin Elmer / Hitachi	200	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	220	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	300	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	320	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	329	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	330	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	340	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	556	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	557	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	220A	PE Pt. 239-0354	D005
Perkin Elmer / Hitachi	330 LC	PE Pt. 239-0354	D005
Perkin Elmer	272	PE Pt. 0057-0194	D013
Perkin Elmer	280	PE Pt. 0057-0194	D013
Perkin Elmer	360	PE Pt. 0057-0194	D013
Perkin Elmer	370	PE Pt. 0057-0194	D013
Perkin Elmer	372	PE Pt. 0057-0194	D013
Perkin Elmer	373	PE Pt. 0057-0194	D013
Perkin Elmer	380	PE Pt. 0057-0194	D013
Perkin Elmer	460	PE Pt. 0057-0194	D013
Perkin Elmer	560	PE Pt. 0057-0194	D013

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Perkin Elmer	575	PE Pt. 0057-0194	D002
Perkin Elmer	603	PE Pt. 0057-0194	D013
Perkin Elmer	703	PE Pt. 0057-0194	D013
Perkin Elmer	2280	PE Pt. 0057-0194	D013
Perkin Elmer	2380	PE Pt. 0057-0194	D013
Perkin Elmer	3030	PE Pt. 0057-0194	D013
Perkin Elmer	3030B	PE Pt. 0057-0194	D013
Perkin Elmer	4000	PE Pt. 0057-0194	D013
Perkin Elmer	5000	PE Pt. 0057-0194	D013
Perkin Elmer	5100	PE Pt. 0057-0194	D013
Perkin Elmer	5100PC	PE Pt. 0057-0194	D013
Perkin Elmer	5100Z Leeman	PE Pt. 0057-0194	D013
Perkin Elmer	EZ 201	PE Pt. C055-5050	D002
Perkin Elmer	EZ 301	PE Pt. C055-5050	D002
Perkin Elmer	EZZ201	Prealigned	D005
Perkin Elmer	Integral 2000	PE Pt. C055-5050	D002
Perkin Elmer	Integral 4000	PE Pt. C055-5050	D002
Perkin Elmer	Lambda 1	PE Pt. C055-5050	D002
Perkin Elmer	Lambda 3	PE Pt. C055-5050	D002
Perkin Elmer	Lambda 3A	PE Pt. C055-5050	D002
Perkin Elmer	Lambda 3B	PE Pt. C055-5050	D002
Perkin Elmer	Lambda 6	PE Pt. C055-5050	D002
Perkin Elmer	LC 55	PE Pt. 087-1171	D002
Perkin Elmer	LC 55B		D002
Perkin Elmer	LC 65B		D002
Perkin Elmer	LC 65T		D002
Perkin Elmer	LC 75		D002
Perkin Elmer	LC 85		D002
Perkin Elmer	LC 85B		D002
Perkin Elmer	LC 95	PE Pt. 0271-1709	D002
Perkin Elmer	LC 135		D002
Perkin Elmer	LC 141 MC		D002
Perkin Elmer	LC 241 MC		D002
Pharmacia	TDS from S/N 206	Resolder old circuit board	D002
Pharmacia	Ultrospec II from S/N 6598	Resolder old circuit board	D002
Pharmacia	Ultrospec III	Resolder old circuit board	D002
Pharmacia	Ultrospec Plus from S/N 176	Resolder old circuit board	D002
Shimadzu	6500AA		D001
Shimadzu	AA-670		D003
Shimadzu	AA-670G		D003
Shimadzu	AA-6200		D003
Shimadzu	AA-6300		D003
Shimadzu	AA-6500		D003
Shimadzu	AA-6800		D003
Shimadzu	CS-920		D003
Shimadzu	CS-930		D003
Shimadzu	CS-9000		D003
Shimadzu	D300L		D003

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Shimadzu	LC 3A		D001
Shimadzu	LC 4A		D001
Shimadzu	LC 6A		D001
Shimadzu	MPS 2000		D003
Shimadzu	SP4		D001
Shimadzu	SPD-2A		D001
Shimadzu	SPD-6A		D001
Shimadzu	SPD-6AV		D001
Shimadzu	SPD-7A		D001
Shimadzu	SPD-7AV		D001
Shimadzu	UV 106A		D003
Shimadzu	UV 120		D003
Shimadzu	UV 120-1		D003
Shimadzu	UV 120-2		D003
Shimadzu	UV 140		D003
Shimadzu	UV 150		D003
Shimadzu	UV 160		D003
Shimadzu	UV 160A		D003
Shimadzu	UV 180		D003
Shimadzu	UV 190		D003
Shimadzu	UV 200		D003
Shimadzu	UV 202		D003
Shimadzu	UV 210		D003
Shimadzu	UV 210A		D003
Shimadzu	UV 240		D003
Shimadzu	UV 250		D003
Shimadzu	UV 260		D003
Shimadzu	UV 265		D003
Shimadzu	UV 350		D003
Shimadzu	UV 360		D003
Shimadzu	UV 365		D003
Shimadzu	UV 370		D003
Shimadzu	UV 625		D003
Shimadzu	UV 1200		D003
Shimadzu	UV 1201		D003
Shimadzu	UV 1202		D003
Shimadzu	UV 1600		D003
Shimadzu	UV 1700		D003
Shimadzu	UV 1800		D003
Shimadzu	UV 2100		D003
Shimadzu	UV 2401 PC		D003
Shimadzu	UV 2501 PC		D003
Shimadzu	UV 2600		D003
Shimadzu	UV 2700		D003
Shimadzu	UV 3000		D003
Shimadzu	UV 3100		D003
Shimadzu	UV 3600		D003
Shimadzu	UV 3700		D003
Shimadzu	UV mini 1240		D003

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia



## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Shimadzu	UV -Vis 160U		D003
Shimadzu	UV-Vis 190		D003
Spectra Physics	3100		D016
Thermo Separation	Spectro Monitor 3200	Use clip from old lamp	D017
Thermo Spectronic	Aquamate	Prealigned	D006
Thermo Spectronic	Evolution 600	Prealigned	D006
Thermo Spectronic	UV1	Prealigned	D006
Thermo Spectronic	UV1-100	Prealigned	D006
Thermo Spectronic	UV2	Prealigned	D006
Thermo Spectronic	UV3	Prealigned	D006
Thermo Spectronic	UV4	Prealigned	D006
Tosoh	8010		D016
Tosoh	8011		D016
Unicam	Helios Alpha	Prealigned	D006
Unicam	Helios Beta	Prealigned	D006
Unicam	Helios Delta	Prealigned	D006
Unicam	Helios Epsilon	Prealigned	D006
Unicam	Helios Gamma	Prealigned	D006
Unicam	Helios II	Prealigned	D006
Unicam	Solar S4	Prealigned	D010
Unicam	Solaar 929		D014
Unicam	Solaar 939		D014
Unicam	Solaar 959		D014
Unicam	Solaar 969		D014
Unicam	SP7500		D017
Unicam	UV1	Prealigned	D006
Unicam	UV2	Prealigned	D006
Unicam	UV2-200	Prealigned	D006
Unicam	UV3	Prealigned	D006
Unicam	UV4	Prealigned	D006
Unicam	UV4-100	Prealigned	D006
Unicam	UV4-200	Prealigned	D006
Unicam	UV4-500	Prealigned	D006
Unicam	UV300	Prealigned	D006
Unicam	UV500	Prealigned	D006
Unicam	UV550	Prealigned	D006
Varian	75 Series	Varian Pt. 56-100218-00	D009
Varian	AA200	Varian Pt. 56-100218-00	D009
Varian	AA275	Varian Pt. 56-100218-00	D009
Varian	AA300	Varian Pt. 56-100218-00	D009
Varian	AA400	Varian Pt. 56-100218-00	D009
Varian	AA475	Varian Pt. 56-100218-00	D009
Varian	AA600	Varian Pt. 56-100218-00	D009
Varian	AA640	Varian Pt. 56-100218-00	D009
Varian	AA775	Varian Pt. 56-100218-00	D009
Varian	AA800	Varian Pt. 56-100218-00	D009
Varian	AA875	Varian Pt. 56-100218-00	D009
Varian	AA975	Varian Pt. 56-100218-00	D009

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia

## Deuterium Lamp Compatibility List

Manufacturer	Model	Comments	Part Number
Varian	AA1275	Varian Pt. 56-100218-00	D009
Varian	AA1475	Varian Pt. 56-100218-00	D009
Varian	AA3000	Varian Pt. 56-100218-00	D009
Varian	Cary 50	Varian Pt. 56-100218-00	D009
Varian	Cary 100	Varian Pt. 56-100218-00	D009
Varian	DMS 100	Varian Pt. 56-100218-00	D009
Varian	LC 9050	Use bracket from old lamp	D018
Varian	ProStar 310	Use bracket from old lamp	D018
Varian	Spectra 800	Varian Pt. 56-100218-00	D009
Varian	Spectra A10	Varian Pt. 56-100218-00	D009
Varian	Spectra A20	Varian Pt. 56-100218-00	D009
Varian	Spectra A30	Varian Pt. 56-100218-00	D009
Varian	Spectra A40	Varian Pt. 56-100218-00	D009
Varian	Spectra A640	Varian Pt. 56-100218-00	D009
Varian	Spectro AA 220FS	Varian Pt. 56-100218-00	D009
Varian	UV Superscan	Varian Pt. 56-100218-00	D009
Varian	UV-50	Use bracket from old lamp	D018
Varian	UV-100	Use bracket from old lamp	D018
Varian	UV-200	Use bracket from old lamp	D018
Varian	UV 634	Varian Pt. 56-100132-00	D012
Varian	UV 635	Varian Pt. 56-100132-00	D012
Varian	Vista 5000	Use bracket from old lamp	D018
Varian	Vista 5500	Use bracket from old lamp	D018
Waters	480LC		D016
Waters	481	Use bracket from old lamp	D016
Waters	481LC		D016
Waters	Lambda max		D016

[www.spectrolamps.com.au](http://www.spectrolamps.com.au)

[sales@spectrolamps.com.au](mailto:sales@spectrolamps.com.au)

PO Box 665, Endeavour Hills, 3802, VIC, Australia